

Instructions to use Zoll Data Control

- Gather patient name and run number
- Insert card into the Omni Drive till the red light flashes on. Do not force the card into the drive. The side of the card with many holes goes into the drive first.
- Double click Zoll Data Control Icon on the desktop
- Click File
 - Read card
- Pop up appears
 - Click Read Card again
- Pop up appears
 - Show percentage of card as it is read
- You may receive an error message
 - Gives error code of zero
 - Click ok till it goes away, may require many clicks
- Goes to processing card with progress bar
- Goes to New Case Summary
 - List all runs downloaded on card
 - Make note of date, time and unit serial number
- Click Create cases
- When finished creating cases
 - The date you need should appear in the left pane
 - If it does not:
 - Click on down arrow beside filter
 - Select appropriate time frame
- Click + sign beside date you need
- Click + beside "Unknown"
- Look for appropriate date, time and serial number
- Click on it
- Case opens
- Select General tab at top of case
 - Enter run number
 - Enter Patient Name
- Do this for each downloaded case
- Click File
 - Erase card
- Pop up appears
 - Click erase card again
- Pop up shows progress
- Remove card by pushing on eject button
- Return card to monitor IMMEDIATELY!!!!!!